

## Certificate of Analysis

Page 1 of 2

<b>Client:</b>	OM Natural Health Limited	<b>Lab No:</b>	4026283	POSSSEP-4v1
<b>Contact:</b>	Ian Ord C/- OM Natural Health Limited 1 Anzac Avenue Whakatane 3120	<b>Date Received:</b>	07-Nov-2025	
		<b>Date Reported:</b>	12-Nov-2025	
		<b>Quote No:</b>		
		<b>Order No:</b>		
		<b>Client Reference:</b>		
		<b>Submitted By:</b>	Ian Ord	

### Sample Type: Plant Derived Food Additives and Supplements

<b>Sample Name:</b>	Reishi Capsules
<b>Lab Number:</b>	4026283.4
Multiresidue Analysis 2 - Type D Samples	
<b>Analytes Detected:</b>	None

Please refer to the detection limits table for the list of analytes screened and their detection limits.

## Summary of Methods

The following table(s) gives a brief description of the methods used to conduct the analyses for this job. The detection limits given below are those attainable in a relatively simple matrix. Detection limits may be higher for individual samples should insufficient sample be available, or if the matrix requires that dilutions be performed during analysis. A detection limit range indicates the lowest and highest detection limits in the associated suite of analytes. A full listing of compounds and detection limits are available from the laboratory upon request. Unless otherwise indicated, analyses were performed at Hill Labs, 28 Duke Street, Frankton, Hamilton 3204.

### Sample Type: Plant Derived Food Additives and Supplements

Test	Method Description	Default Detection Limit	Sample No
Multiresidue Analysis 2 - Type D Samples	Solvent extraction, SPE cleanup, dilution. Analysis by LC-MS/MS. In-house (using a Citrate buffered QuEChERS extraction).	0.10 - 2 mg/kg	4

These samples were collected by yourselves (or your agent) and analysed as received at the laboratory.

Testing was completed on 12-Nov-2025. For completion dates of individual analyses please contact the laboratory.

Samples are held at the laboratory after reporting for a length of time based on the stability of the samples and analytes being tested (considering any preservation used), and the storage space available. Once the storage period is completed, the samples are discarded unless otherwise agreed with the customer. Extended storage times may incur additional charges.

This certificate of analysis must not be reproduced, except in full, without the written consent of the signatory.



Shaun Clay BSc  
Senior Technologist - Food and Bioanalytical

## Detection Limits

Analytes	Detection Limit	Analytes	Detection Limit	Analytes	Detection Limit
Multiresidue Analysis 2 - Type D Samples		Flusulfamide	0.10 mg/kg	Spirotetramat-cis-keto-hydroxy	0.10 mg/kg
Abamectin	0.10 mg/kg	Flutolanil	0.10 mg/kg	Spirotetramat-enol-glucoside	0.10 mg/kg
Acetamidiprid	0.10 mg/kg	Fluxapyroxad	0.10 mg/kg	Spirotetramat-mono-hydroxy	0.10 mg/kg
Aldicarb	0.10 mg/kg	Forchlorfenuron	0.10 mg/kg	Spiroxamine	0.10 mg/kg
Aldicarb sulfone	0.10 mg/kg	Gibberellic acid (GA3)	0.10 mg/kg	Sulfoxaflor	0.10 mg/kg
Aldicarb sulfoxide	0.10 mg/kg	Halauxifen-methyl	0.10 mg/kg	Tebufenozide (Mimic)	0.10 mg/kg
Ametoctradin	0.10 mg/kg	Heptenophos	0.10 mg/kg	Teflubenzuron	0.2 mg/kg
Ametryn	0.10 mg/kg	Imidacloprid	0.10 mg/kg	Tepraloxymid	0.10 mg/kg
Anilazine	1.0 mg/kg	Indaziflam	0.10 mg/kg	Tetraconazole	0.10 mg/kg
Anilofos	0.10 mg/kg	Iodocarb (IPBC)	0.2 mg/kg	Thiabendazole	0.10 mg/kg
Azadirachtin	0.2 mg/kg	Ipconazole	0.10 mg/kg	Thiacloprid	0.10 mg/kg
Benzalkonium chloride (C10)	0.3 mg/kg	Isoproturon	0.10 mg/kg	Thiamethoxam	0.10 mg/kg
Benzalkonium chloride (C12)	1.0 mg/kg	Isopyrazam	0.10 mg/kg	Thifluzamide	0.10 mg/kg
Benzalkonium chloride (C14)	1.0 mg/kg	Isoxathion	0.10 mg/kg	Thiophanate-methyl	0.10 mg/kg
Benzalkonium chloride (C16)	0.3 mg/kg	Isoxathion oxon	0.10 mg/kg	Triadimenol	0.15 mg/kg
Bixafen	0.10 mg/kg	Lufenuron	0.7 mg/kg	Tribenuron-methyl	0.10 mg/kg
Boscalid	0.10 mg/kg	Mandipropamid	0.10 mg/kg	Trichlorfon	0.10 mg/kg
Carbendazim (including Benomyl and Thiophanate)	0.10 mg/kg	Mandestrobin	0.10 mg/kg	Triflumuron	0.10 mg/kg
Carfentrazone-ethyl	0.10 mg/kg	Mefentrifluconazole	0.10 mg/kg	Triforine	0.10 mg/kg
Chlorantraniliprole	0.10 mg/kg	Metamitron	0.10 mg/kg	Uniconazole	0.10 mg/kg
Chloridazon	0.10 mg/kg	Metaldehyde	0.3 mg/kg		
Clethodim	0.10 mg/kg	Metconazole	0.10 mg/kg		
Clofentezine	0.2 mg/kg	Methabenzthiazuron	0.10 mg/kg		
Clothianidin	0.10 mg/kg	Methomyl	0.10 mg/kg		
Cyantraniliprole	0.10 mg/kg	Methoxyfenozide	0.10 mg/kg		
Cyazofamid	0.10 mg/kg	Methyl anthranilate	0.2 mg/kg		
Cyflufenamid	0.10 mg/kg	Metrafenone	0.10 mg/kg		
Cymoxanil	0.10 mg/kg	Milbemectin	0.2 mg/kg		
Didecyldimethylammonium chloride (DDAC)	0.10 mg/kg	Nicosulfuron	0.10 mg/kg		
Desmedipham	0.10 mg/kg	Novaluron	0.10 mg/kg		
Diethofencarb	0.10 mg/kg	Octhilinone	0.2 mg/kg		
Difflubenzuron	0.2 mg/kg	Oryzalin	0.10 mg/kg		
Dinotefuran	0.10 mg/kg	Oxamyl	0.10 mg/kg		
Dodine	0.7 mg/kg	Oxathiapirolin	0.10 mg/kg		
Emamectin	0.3 mg/kg	Pencycuron	0.2 mg/kg		
Empenthrin	2 mg/kg	Penflufen	0.10 mg/kg		
Ethofumesate	0.10 mg/kg	Penthiopyrad	0.10 mg/kg		
Etobenzanid	0.10 mg/kg	Phenmedipham	0.10 mg/kg		
Fenamidone	0.10 mg/kg	Propamocarb	0.10 mg/kg		
Fenbuconazole	0.10 mg/kg	Propargite	0.10 mg/kg		
Fenhexamid	0.3 mg/kg	Proquinazid	0.10 mg/kg		
Fenoxycarb	0.10 mg/kg	Prosulfocarb	0.10 mg/kg		
Fenpropidin	0.10 mg/kg	Prothioconazole-desthio	0.10 mg/kg		
Fenpyrazamine	0.10 mg/kg	Pydiflumetofen	0.10 mg/kg		
Fenpyroximate	0.10 mg/kg	Pyraclostrobin	0.10 mg/kg		
Fipronil	0.10 mg/kg	Pyridaphenthion	0.10 mg/kg		
Fonicamid	0.10 mg/kg	Pyrifluquinazon	0.10 mg/kg		
Fluazinam	0.10 mg/kg	Pyriofenone	0.10 mg/kg		
Flufenacet	0.10 mg/kg	Quinoxifen	0.10 mg/kg		
Flufenoxuron	0.10 mg/kg	Saflufenacil	0.10 mg/kg		
Florasulam	0.10 mg/kg	Sethoxydim	0.10 mg/kg		
Flumethrin	0.10 mg/kg	Spinetoram	0.10 mg/kg		
Flumetsulam	0.10 mg/kg	Spinosad	0.10 mg/kg		
Flumioxazin	0.10 mg/kg	Spiromesifen	0.10 mg/kg		
Fluopicolide	0.10 mg/kg	Spiromesifen-enol	0.10 mg/kg		
Fluopyram	0.10 mg/kg	Spirotetramat	0.10 mg/kg		
		Spirotetramat-cis-enol	0.10 mg/kg		